

# QUICKRAY<sup>PRO</sup>



## NEW GENERATION Intraoral Dental X-RAY SENSOR

Unmatched Image Quality, Unbelievable Value

---

The **Quickray<sup>PRO</sup>** digital x-ray sensor offers the best resolution on the market, backed by a market leading warranty. For high precision diagnosis, patient comfort & reliability choose **Quickray<sup>PRO</sup>**!



TWAIN compatible, Quickray<sup>PRO</sup> works in all leading imaging software and can be used in conjunction with existing sensors.

COMPATIBLE WITH






# Your Practice Deserves To Purchase The Right Sensor! Let Our Team of Experts Help You Make the Best Choice

## SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs=35.7 lp/mm
SENSOR SIZES	Size 1: 36.2 x 24.1 x 5.9 mm Size 2: 42.8 x 30.5 x 5.7 mm
ACTIVE AREA	Size 1: 29.5 x 19.7 mm Size 2: 36.0 x 25.8 mm
PIXEL SIZE	14x14um
COMPATIBILITY	Windows 7/ 10, USB 2.0
BODY THICKNESS	5.7MM
SCINTILLATOR	CSI



## QUICKRAY<sub>PRO</sub> SENSOR FEATURES



ACTUAL IMAGE



### CORD

- Optimized image quality at low dose
- Comfortable design
- Ease of integration
- Outstanding reliability
- Available in size 1 & 2



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity
- 100,000+ bendings RA test passed

## SENSOR KIT INCLUDES



Protective Sheaths



Uni-Verse-All Positioning Kit



Wall Mount Holder



## Warranty

5 Year Warranty  
7 Year Replacement Plan



## Expedited Rental Service

Available for NO Downtime!

**CONTACT US TODAY!**



Outstanding Dental Equipment Since 1989

www.VideoDental.com | sales@videodental.com | 800.323.2690

156 East Granada Boulevard, Ormond Beach, FL 32176