

iCAM



The **Low Cost, High Quality** Intraoral Camera

Perfect for Hygiene and Patient Education

The camera is easy to use and works with most imaging software. Unlike most inexpensive cameras which have a short life span, render a mediocre image, and whose capture button does not save an image in your software; the iCam performs admirably, delivering an excellent image, with its programmable capture button it is compatible with all major dental imaging software.



iCam works in all leading imaging software and can be used in conjunction with existing cameras.

COMPATIBLE WITH



Your Practice Deserves To Purchase The *Right* Camera!
 Let Our Team of Experts Help You Make the Best Choice

TECHNICAL SPECIFICATIONS

ANGLE OF LENS	90 Degrees
CAMERA FOCUS	Fixed
DEPTH OF FIELD	40mm
LIGHTING SOURCE	6 Lights Emitting Diodes
WARRANTY	1 Year Parts and Labor
FOOTPRINT	Less Than 6 sq Inch
INTERFACE	Universal Serial Bus (USB)
RESOLUTION	1.3MP



CAMERA KIT INCLUDES



Protective Sheaths



iCam



Wall Mount Holder



Warranty

1 year parts and labor

**CONTACT YOUR
DEALER TODAY!**



© Copyright 2018 DENTERPRISE INTERNATIONAL, INC. All Rights Reserved. Information Subject to Change Without Notice.

[DEALER TAG HERE]

Denterpriseintl.com