

YOUR HANDHELD HERO

ECOTRON DT-703



***Be an imaging HERO with
our handheld marvel.***

ECOTRON... YOUR HANDHELD HERO

***The most advanced
lightweight x-ray device
available on the market!***

- ✓70kV power for better image quality and cavity detection.
- ✓3mA for quick exposure.
- ✓0.3mm focal spot for sharper images.
- ✓3-second exposure intervals that allow faster FMX and BW series.
- ✓Device design allows single-handed ambidextrous operation
- ✓Lightweight of 4.7 lbs. for better ergonomics and to prevent user fatigue.
- ✓The soft shoulder strap for drop prevention, and improved comfort.
- ✓White glove warranty service: 2-year warranty and free loaners.



Visit Our Website

www.globalimagingusa.com



Handheld x-rays are easy to handle and faster to capture FMX. Rapid duty cycle of 3 seconds guarantees the x-ray will always be ready for your next exposure. Patient movement can be observed and image acquisition is quicker without leaving the room. Soft neck strap allows efficient sensor repositioning and prevents drops. Freehanded aiming is easier than positioning wall mounted x-ray, especially when traditional x-ray arm is drifting.

Independent testing found that handheld x-rays are safe and easy to use, and radiation levels are well below FDA limits. Image quality is comparable to wall-mounted x-rays. Always practice ALARA principle when imaging.

sales@globalimagingusa.com

(702) 850-2797