


Extraoral Dental Suction Unit

The Extraoral Dental Suction Unit provides up to 99.9% air purification

With 4 filtration levels and UV sterilization, the Extraoral Dental Suction Unit kills 99.9% of coronaviruses, providing increased comfort and peace of mind for your patients and staff.

- 4 filtration levels, high efficiency closed purification
- Kills 99.9% of coronaviruses
- Effectively removes aerosols and particulate matter
- In many cases, allows you to work without an assistant
- Great for hygiene use
- High-speed, super-quiet performance
- 360-degree multi-angle operation for ease of use
- Durable construction

Four-level high-efficiency closed purification



Primary Filter: Separates the inhaled water mist and filters larger particles, dust, and more.

Sterilization Filter: Filters microorganisms and small particles.

UV Plasma Sterilization: Disperses ultraviolet plasma organic molecules to effectively kill viruses and bacteria.

H14 HEPA Filter: Isolates 0.01 micron particle matter.

Flexible suction arm can be adjusted to any angle with optimum stability

