

Product Manual

Product Name 产品名称	Rabbit Anti-Goat IgG(H+L) Polyclonal Antibody, FITC Conjugated 兔抗山羊IgG(H+L)多克隆抗体 -FITC
Catalog No./货号	FK-0166
Buffer/缓冲液	50% glycerol, 0.01 M PBS, 0.5% RSA, 0.1%NaN ₃
Storage/存储	Shipped at 2-8°C. Upon delivery aliquot and store at -20±5°C. Avoid freeze/thaw cycles. Away from light. 在2-8°C运输。收货后，在-20±5°C避光保存，避免反复冻融。
Usage Recommendations 使用建议	Flow-Cyt/ IF:1:150-1500, IHC is recommended at 1:500-2500 and this product can also be used for ELISA. Flow-Cyt/ IF:1:150-1500, IHC建议在1:500-2500，本产品也可以用于ELISA。
Expiry Date/有效期	Please ensure that the product is used within the expiration date indicated on the tube. 请在管体注明有效期内使用。
Matters Needing Attention 注意事项	This product contains preservative. Standard Laboratory Practices should be followed when handling this material. 本品含有防腐剂。处理该产品时应遵循标准实验室操作规程。
Note/备注	Centrifuge before opening to ensure vial contents are completely collected. 开启前请离心，确保完全收集瓶内的产品。

support@fantibody.com

FOR RESEARCH USE OR FURTHER MANUFACTURING ONLY