

# Bestneb Nebuliser Pump

## 12LPM Ultra Quiet Operation

### High quality aerosol compressor

BESTNEB Nebuliser is a high quality aerosol compressor with contemporary, non-medical styling.

It is appealing to the clinician and the patient. Especially in the pediatric case, it makes children's respiratory therapy less frightening and has revolutionised the industry's approach to children's home health care products. According to the latest study, the majority of nebuliser user are from the paediatric field due to the poor resistance to air pollution.

### Features

- Contemporary styling
- Ultra quiet operation
- Pleasant pastel colours
- Compartment for storing accessories
- Made from high quality polymousse
- High quality aerosol compressor
- Large air intake filter



## Product Information

### Specifications

Power: 220/240 V ~ 50 HZ

Dimensions: 29 x 20 x 12 cm

Weight: 2.5 Kg

Maximum pressure: 32 PSI

Flow: 12 LPM

Sound level: 50 DB

Particle size (mmad): 1.1 micron

Warranty: 2 Years



Code	Description
PNB	Bestneb nebuliser pump
PNLF	Spare filter (PK5)
31310	Nebuliser tubing 2 metres
3569	Nebeasy Nebuliser bowl, mask and tubing (set) - Adult
3570	Nebeasy Nebuliser bowl, mask and tubing (set) - Child