


# A simple, yet effective solution for a continuous spray system.



## OtoClear Tips | Key Features

- 

**Gentle Touch Design™**  
Improved patient comfort with soft, gentle design.
- 

**Three Stream Design**  
Directs water flow away from the tympanic membrane.
- 

**Exit Portals**  
Prohibits pressure build-up and allows for direct drainage into an ear basin, reducing cleanup time.
- 

**Flared Tip Design**  
Eliminates over-insertion into the ear canal.

## OtoClear®, powered by Aquabot®.

Aquabot is a simple, yet effective pressure spray bottle system that attaches the OtoClear Ear Irrigation Tips at the end of the device wand creating a continuous and safe spray for ear irrigation procedures. Healthcare professionals can operate this system with ease through a simple pump and push release motion that offers an uninterrupted spray process.

### Otoclear Product Offerings Overview

#### Aquabot Reorder items



**#7200**  
OtoClear Tips  
OtoClear tips come with 40 single-use tips per box.



**#3700**  
Ear Basins  
Bionix Ear Irrigation Basins are sold in a box of three and available to be purchased individually.



\*To reorder OtoClear Tips, see product offering #7200

MKT-0077  
Rev.A

#### Aquabot Product Offerings

##### #7270

- Aquabot
- 20 Tips
- 1 Ear Basin

##### #7275

- Aquabot
- 5 Tips
- 1 Ear Basin

MEDICAL PRODUCTS



WESTERN AUSTRALIA

14 Stirling Highway, Nedlands  
Western Australia, 6009

 (08) 9273 3022

 (08) 9273 3017

 sales@amawa.com.au