




BLANKET WARMERS

NULINE
REFRIGERATION

FEATURES

- Fan Forced System
- Overheating Alarm
- Safety Cutoff at 70°C
- Double Glazed Doors
- Stainless int./ext. with Chrome Plated Shelves
- Microprocessor Control
- Digital Display with High Alarm
- TGA Certified



	60°C
	80 L 247 L 611 L 1000 L
	500w x 600d x 625h 700w x 700d x 920h 750w x 700d x 1880h 1200w x 700d x 1880h

+61 3 9562 3355
info@nuline.au
7 Longford Court, Springvale VIC 3171
www.nulinerefrigeration.com.au



BLANKET WARMERS



Model	BW100	BW260	BW650	BW1100
Gross Capacity	80 L	247 L	611 L	1000 L
Blanket Capacity	5 - 6	10 - 12	30 - 35	47 - 52
External Dimensions in mm (w/d/h including castors)	500 x 600 x 625	700 x 700 x 920	750 x 700 x 1880	1200 x 700 x 1880
Internal Dimensions in mm (w/d/h)	390 x 480 x 370	590 x 685 x 600	630 x 610 x 1590	1100 x 550 x 1555
Packed Dimensions in mm (w/d/h including pallet)	575 x 670 x 765	800 x 800 x 1070	850 x 800 x 2030	1300 x 800 x 2030
Gross Weight - kg (packed weight)	63	109	180	260
Net Weight - kg	48	89	160	220
Voltage	240V			
Temperature Range	60°C			
Exterior Cabinet Finish	Stainless Steel			
Inner Liner Material	Stainless Steel			
Door	Glass			
Door Hinge		RH		LH / RH
Door Handle	Full Length			
Type of Controller	Microprocessor			
Temperature Display	Digital			
BMS Connectivity				
(BMS Connectivity for Remote Alarms only)	Optional			
(BMS Connectivity for 4-20mA Transmitter)	Optional			
Port Entry	Optional			
Power Requirement	10amp GPO			
Castors	Fixed Legs		Yes	
Adjustable Shelves	2	3		6



Nuline Refrigeration reserves the right to alter the specifications without prior notice.

BLANKET WARMERS

NULINE
REFRIGERATION

Since 1985, Nuline Refrigeration has been manufacturing refrigerators and freezer cabinets for the Australian market. Our specialised range of biomedical cabinets are used extensively across the hospitals, medical centres, pharmacies and laboratories.

Focusing on customer needs, we manufacture, import and customise the biomedical cabinets from a range of selected global suppliers. Our factory is based in Melbourne, Australia and we are proudly 100% Australian owned business.

Each of our manufactured cabinets are custom built based on customer orders by our craftsmen right here in Melbourne.

We have extensive capability to design, engineer and build a range of refrigerated cabinets for you. Our standard product range includes, fridges, freezers (up to -40°C), blanket warmer, mortuary refrigerators and refrigerated incubators.

