



## DIAGNOSTIC LITHIUM ION MACROVIEW SET

STOCK CODE: WEL97206MVS

- **Macroview Otoscope with throat Illuminator Head**
  - A nearly complete view of the tympanic membrane with approximately twice the field of view and 30% greater magnification than a traditional otoscope
  - Better clarity & definition of landmarks
  - The ability to adjust focus for variable ear canal length or farsighted eyes
  - Tip grip for secure fastening and ejection of ear specula
  - Fiber optics produce cool light with no reflections, no obstruction
- **Coaxial Ophthalmoscope Head**
  - Halogen HPX™ lamp provides 30% whiter, brighter light than standard halogen lamps for true tissue colour & long-lasting performance
  - Coaxial optics produce a shadow-free spot, easier entry into undilated pupils, and a larger field of view versus standard ophthalmoscopes
  - Detect corneal abrasions with cobalt blue filter
  - Polarizing filter virtually eliminates corneal reflection
  - Sealed optics keep out dust & dirt
  - 18 unique aperture/filter combinations for greater versatility
  - Red-free filter may be used with any aperture
  - 28 focusing lenses with a range of -25 to +40 diopters
- Lithium Ion Rechargeable Handle, Pod & Charger
- Twice the on-time as a conventional handle, half the weight & fits more comfortably in your hand
- Black carry case