

# 3-L Precision Syringe

For peace of mind regarding the accuracy of your spirometer

The 3-L Precision Syringe from Vitalograph® allows quick, easy accuracy checks for all spirometers, as recommended by the ATS/ERS Standardization of Spirometry 2019 Update and required by the Spirometry standard ISO 26782:2009.

The 3-L Precision Syringe provides reliable verification of accuracy and can be used with both flow or volume sensor systems and is compatible with any manufacturers' equipment. An in-line filter adapter is included to allow use of Bacterial-Viral Filters for protection against contamination.

## Technical Specifications

|                               |                         |
|-------------------------------|-------------------------|
| <b>Model Number:</b>          | 2040/36020              |
| <b>Volume:</b>                | 3L                      |
| <b>Accuracy:</b>              | 0.5%                    |
| <b>Operating Temperature:</b> | 15–40°C                 |
| <b>QA/GMP Standards:</b>      | ISO 13485, FDA 21CFR820 |
| <b>Gross Size:</b>            | 500mm x 170mm x 125 mm  |
| <b>Gross Weight:</b>          | 2.5kg                   |



Vitalograph®, Vitalograph® are registered trademarks of Vitalograph Ltd.



Distributed by:  
AMA Medical Products  
amamedicalproducts.com.au  
sales@amamedicalproducts.com.au  
1800 626 292

A Global Leader in Respiratory Diagnostics