



Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Classic II Pediatric and Infant Stethoscopes	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3200
	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.			Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more of what you need to hear with amplification of 24x, ambient noise reduction and optional connectivity.
Approximate Length	28"/71cm	27"/68cm	27"/69cm	22"/56cm 27"/69cm	22"/56cm 27"/68cm	27"/68cm
Approximate Weight	105g, 95g	160g	150g	167g, 177g	175g, 185g	185g
Warranty	3 year	3 year	5 year	7 year	7 year	2 year
Diaphragm	Non-Tunable (Floating)	Tunable	Tunable (2)	Tunable (2)	Tunable	Electronic
Approximate Chestpiece Diameter	Double-sided Pediatric 1.4"/3.4cm 1.1"/2.7cm Infant 1.1"/2.8cm 0.7"/1.7cm	Single-sided 1.7"/4.3cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.9"/4.8cm	Single-sided 2"/5.1cm
Features	Specially designed for the smallest patients	Designed for easier grip and maneuvering	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Eliminates 85% (on average) of ambient noise Bluetooth® wireless technology Recording and playback



The power to learn, in the palm of your hand.

3M™ Littmann® Learning Institute App – Free with purchase of a 3M™ Littmann® Stethoscope!



3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

Components
(Chestpiece/Eartube/Stem)

Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Model 3200 Stethoscope
-------------------------------	----------------------------------	--------------------------------	-------------------------------------	---------------------------------------	--------------------------------	-----------------------------------

All Standard Finish	2114 Black	2113 Black	2144L Black	5620 Black	5633 Pearl Pink	6152 Black	6156 Plum	2160 Black	3200BK27 Black
	2214R Red	2113R Red	2146 Burgundy	5621 Grey	5831 Plum	6154 Navy	6158 Raspberry	2163 Burgundy	3200BU Burgundy
	2124 Caribbean	2119 Caribbean	2147 Navy	5623 Caribbean	5832 Lavender	6155 Hunter Green		2164 Navy	3200NB Navy
		2122 Raspberry		5627 Burgundy	5835 Turquoise			2167 Plum	
			5622 Navy	5839 Lemon-Lime					
			5630 Ceil Blue						
All Black			2141 Black	5803 Black	5868 Burgundy	6163 Black		2161 Black	
						6168 Navy			
All Brass Finish									
All Copper Finish				5809 Chocolate					
All Smoke				5811 Black	5871 Raspberry	6162 Black		2176 Black	
				5812 Dark Olive		6166 Plum			
Rainbow/Brass/Brass				5806 Raspberry	5807 Caribbean		6171 Turquoise		
Brass/Smoke/Smoke		2153 Caribbean				6164 Black			
Rainbow/Black/Black				5870 Black		6165 Black			
Mirror/Smoke/Smoke				5862 Raspberry	5863 Navy				
Mirror/Standard/Mirror						6170 Burgundy	6177 Black		
Champagne/Smoke/Smoke				5861 Black	5864 Burgundy	6179 Black			
Champagne/Black/Champagne						6176 Burgundy	6190 Caribbean		
COLORED STEM SERIES: (square inset represents stem color)									
Black/Black/Color						6200 Black	6203 Black		
						6201 Black			
High-Polish Smoke/Black/Color						6202 Navy	6204 Black		
Rainbow/Black/Color						6205 Plum			
Champagne/Champagne/Color						6206 Hunter			
Smoke/Smoke/Color				5872 Turquoise	5875 Lime				
				5873 Gray					
Mirror/Standard/Color				5874 Caribbean					

Components:
Chestpiece/Eartube/Stem

Standard Finish

Black

Brass Finish

Copper Finish

Smoke

Rainbow

Mirror

Champagne

High-Polish Smoke