

LKUv 1610

Laboratory refrigerator with
Comfort controller
MediLine



Gross/net capacity	142/133 litre
Ext. dimensions in mm (w/d/h)	601/618/820
Int. dimensions in mm (w/d/h)	440/435/670
Energy consumption in 365 days	273 kWh
Ambient temperature range	+10 °C to +35 °C
Refrigerant	R 600a
Connection Rating	1.0 A
Frequency / voltage	50 Hz / 220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	+3 °C to +16 °C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	electronic controller
Temperature display	external digital
Alarm in case of malfunction	optisch und akustisch
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	3
Shelf loading	45 kg
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	40 mm
Gross/net weight	40/37 kg
EAN-Nr.	9005382237694

