

LGPv 1420

MediLine



Gross/net capacity	1366/980 litre
Ext. dimensions in mm (w/d/h)	1430/830/2160
Int. dimensions in mm (w/d/h)	1236/650/1550
Energy consumption in 365 days	2654 kWh
Ambient temperature range	+10°C to +40°C ¹
Refrigerant	R 290
Connection Rating	4.5 A
Voltage	220-240 V-
Cooling system	forced-air
Defrosting	automatic
Temperature range	-9°C to -26°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	304-grade stainless steel
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	8
Shelves, freezer compartment	8
Castors	castors with brake at the front, castors at the rear
Lock	fitted
Self-closing door	yes
Door hinges	left/right
Insulation	83-60 mm
Defrost water drain	yes
Gross/net weight	253/207 kg
EAN-Nr.	9005382197257

¹ Indicates the ambient temperature at which the device can be operated safely.

