




The face of the dial is especially large to make reading the figures much easier.



The rugged metal housing withstands years of use under the hardest conditions imaginable.

## seca 762

Mechanical personal scales

The precise design classic has a large clear and well defined dial. With a large, low-level, non-slip platform and scratch-resistant finish it is easy for the patient to step on and weigh.

- Capacity: 150 kg
- Graduation: 500 g
- Optional: Carrying case seca 421

### Accessories

## seca 421

Carrying case



With this carrying case of water-repellent material, the following seca scales can be safely transported: seca 762 and seca 750.



The large dial makes it easy to read results.

## seca 750

Mechanical floor scales

Pretty and stable: with their pure, light design, these scales will make a good impression in any clinic. A stable, powder-coated steel casing, a circular display and an easy to clean tread area of imitation leather – what more does one need for comfortable weighing? The seca 750 is available in three versions: kilograms, pounds and kilograms/pounds.

- Capacity: 150 kg
- Graduation: 1,000 g
- Optional: Carrying case seca 421

The stable material on the platform ensures safety from slippage.

The housing is of metal treated with powder-coated lacquer.

