

**The Swann-Morton Surgical Blade Remover**

Swann-Morton's Surgical Blade Remover is for single use, and is available sterile in boxes of 50, and non-sterile in boxes of 100. It is designed to remove surgical blades from standard No. 3 and 4 fitment handles, and protects hospital staff removing scalpel blades after operations, guarding against cuts and abrasions.

**Skin Graft Blade (SM9940)**  
Fits Braithwaite, Watson and Cobbett skin graft handles/knives. Box / 10 sterile blades 158mm long.

**'Minor' Disposables**  
Can be used as part of a "Ward Kit". Small plastic handle can be supplied fitted with blade No's 10, 11, 15 or 15A. Polythene blade guard. Non-sterile and bulk-packed.

**Disposable Scalpels**  
Available with all standard stainless steel blade shapes, complete with plastic handles, individually wrapped in sterile packs of 10.

**Guarded Disposable Scalpels**  
With an ultrasonically-welded blade fitment, these scalpels are available with all standard stainless steel blade shapes, and are fitted with a polythene protective guard. Available in a bulk, non-sterile format.

**Disposable Scalpels with Retractable Blade**  
These disposables have been specially designed to reduce accident risk potential by having a high-visibility red handle, retractable blade operated by a colour-coded button which audibly clicks to indicate open and locked position, and ribbed, slip-resistant finger grips for secure hold in gloved hands. They are available with stainless steel blade No's 10 (black button), 11P (blue button) and 15 (green button), in sterile packs of 25. Also available in non-sterile packs of 500.

**Standard Surgical Blades and Handles - No.3 Fitment**  
Surgical Blade No's 6, 9, 10, 10A, 11, E/11, 12, 12D, 14, 15, 15A, 15C, 16, fit No's 3, 3L, 5, 7, 9, B3 & B3L handles. Handles are available in nickel, with selected sizes available in stainless steel. All handles are individually wrapped. The easily gripped B3 & B3L are light, adaptable stainless handles, supplied individually in plastic wallets.

**Standard Surgical Blades and Handles - No.4 Fitment**  
Surgical Blade No's 20, 21, 22, 22A, 23, 24, 25, 25A, 26, 27, fit No's 4, 4L and 6 handles. Handles are available in nickel, with selected sizes available in stainless steel. All handles are individually wrapped.

*Standard surgical blades in stainless steel are sterile and boxed in 100's individually foil-wrapped.*

**'Fine' Blades and Handles**  
The 'Fine' range has small, slim, precise blades in stainless steel for delicate operations. Four special handles with independently screwed top sections to facilitate blade-changing, take all blades, supplied in individual metal foil packets, available in dispenser units of 25. Handles are supplied individually in plastic wallets.

**Swann-Major Surgical Blades and Handles**  
All these extra-strong blades fit either handle. Blades are carbon steel, sterile, individually foil-wrapped, boxed in 50's.

**Post-Mortem Knife - PM40 and PM40B**  
Post-mortem knife with a choice of heavy-duty carbon steel blades. Screw-type fitting to stainless steel handle. Spare blades, individually wrapped, boxed in 10's.

**Stitch Cutters**  
Requiring no handles, the standard stitch cutter is in carbon steel, the midi and long versions in stainless steel. Available individually foil-wrapped in sterile packs of 100.

**BSN medical**  
DISTRIBUTED BY  
Healthcare Division  
Smith & Nephew Pty Ltd  
315 Ferntree Gully Road  
Mount Waverley  
Victoria 3149  
Australia

T: 61 3 8540 6777  
F: 61 3 9544 5086  
www.smith-nephew.com  
Customer Service  
T 13 13 60  
F 1800 671 000

Sterilisation is by Gamma Radiation (minimum dosage 25 K. Gy. 2.5 M. Rad) carried out under strict laboratory control in Swann-Morton (Services) Ltd's Cobalt 60 continuous irradiation plant.

All items are shown actual size

**CE**  
All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC).

**SGS**  
BS EN ISO 9001:1994  
Certificate Number: FM 14713

Re order code: SN 4660

# SURGICAL BLADES AND HANDLES BY

# Swann-Morton®



**Myringotomy Blade**  
Long, angled, stainless steel blade with short cutting edges at tip. Purpose-designed for both up and down cutting. Individually packed, sterilised and boxed in 10's. This blade fits the 'FINE' range of handles.

**Cervical Biopsy Blade**  
In stainless steel, sterilised, supplied in single plastic trays in boxes of 10. Fits both Swann-Major handles.

PM40

PM40B

PM40 (Stainless Steel)