

ALPHA-TCEM™ NE TEMPORARY CEMENT, EUGENOL FREE

ALPHA-TCEM™ NE PROVIDES \downarrow

NON-EUGENOL

STRONG ADHESION

AUTOMIX CARTRIDGE

EASY REMOVAL

EASY TO USE



TEMPORARY CEMENT EUGENOL FREE

Alpha TCEM [™]**NE** is a non-eugenol temporary cement. It is designed for cementing temporary crowns and bridges, inlays and onlays, as well as other temporary restorations. It offers strong adhesion to the tooth yet allows easy removal of the temporary restoration. The quick, automix cartridge makes **Alpha-TCEM** [™]**NE** fast and easy to use.



ALPHA-TCEM™ NE

REF Code 322.0000.002

Product Description 2 x 7g Alpha TCEM™ NE Cartridges, 10 Disposable Mixing Tips

Physical/Mechanical Properties

Working Time Setting Time Compressive strength Film Thickness Acid-soluble arsenic content ≥2 min. at 23 °C ≥5 min. 30 sec. at 37 °C ≤35 Mpa at 24 h ≤25 **µ**m ≤2 mg/kg

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]

1.847.677.5500 PHONE 1.847.677.5502 FAX