



KELLER INTERNATIONAL, LLC


1431 LAKE RD, HAMLIN NY 14464

Info@kellerinternational.com

T: 1-800-578-8772 • F: 1-585-964-8093

KELLER4SALON.COM

UV Hot Towel Warmer – TD207



Function:

This towel warmer is used for heating dry and damp towels, warming cosmetics, storing sterile instruments and accessories (sponges, brushes, napkins, etc.). The sterile environment will hold about two-dozen small hand towels (when rolled up) while maintaining a constant temperature of approximately 176°F.

NOTE: This Towel Warmer is NOT a Sterilization Unit. It will ONLY maintain a sterile environment for ALREADY sterilized equipment.

NOTE: Please follow all your Local / State regulations regarding sterilization. Additional requirements may be necessary.