

# ELECTRIC WHEELCHAIRS



## JRWD503 ECONOMY DUAL FUNCTION POWER WHEELCHAIR

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Dual Functionality allows the wheelchair to operate as an electric wheelchair or as a manual wheelchair.
- Magnesium alloy large rear wheels allow you to reverse safely across obstacles e.g. road humps.
- Aluminum alloy frame makes it very strong, yet lightweight with the wheelchair weighing only 30 kg without batteries.
- Flip-up armrests allow you to get closer to tables and makes it easier to transfer.
- Compact design with foldable seat and backrest, detachable and height-adjustable footrests, allows for easy folding and storage such as the trunk of a car.

## SPECIFICATIONS

Overall Width	: 64 cm (16" seat) or 69 cm (18" seat)
Overall Length	: 108 cm
Overall Height	: 94 cm
Frame Weight	: 30 kg
Battery Weight	: 14 kg
Overall Weight	: 44 kg
Seat Width	: 40 cm or 45 cm (16" or 18")
Seat Depth	: 43 cm (17")
Seat Height	: 49 cm (19")
Max Speed	: 6 km/h
Driving Range	: 15 - 20 km
Rear Wheel Diameter	: 24"
Front Wheel Diameter	: 8"
Turning Radius	: 75cm (30")
Rear Wheel Diameter	: 24"
Front Wheel Diameter	: 8"
Maximum Capacity	: 120kg
Batteries	: 2x12V 20AH SLA
Charger	: 2A c/w Safety Mark
Motors	: 250w