

## **ELECTRIC WHEELCHAIRS**

## ULTRA-LITE AIR ELECTRIC WHEELCHAIR

- LTA compliant personal mobility aid (PMA).
- Each wheelchair comes with 1 quick-release lithium battery. Tools not required.
- Lightweight, only 12.5 kg without battery.
- Each battery weighs 0.8 kg and the chair holds up to 2 batteries at a time.
- 24V, 5.2 AH lithium battery that travel up to 10 km.
- Top speed of 6 km/h.
- Flip-up footrests make it easier for patient to get on or off the chair.



SPECIFICATIONS	
Weight Capacity	: 100 kg
Frame Weight (w/o battery)	: 12.5 kg
Battery Weight	: 0.8 kg x 1 pieces
Total Weight with 1 battery	: 13.3 kg
Total Weight with 2 batteries	: 14.1 kg
Overall Length	: 90 cm
Overall Width	: 67 cm
Overall Height	: 88 cm
Minimum Folded Dimensions	: 66 cm x 35 cm x 77 cm
Seat Size (W x D x H)	: 48 cm x 45 cm x 47 cm
Armrest Height	: 21 cm
Seat-to-Floor Height	: 51 cm
Travelling Range	: 10 km per battery
Max Speed	: 6 km/h
Turning Radius	: 29"
Batteries	: 1 x 24V, 5.2 AH Lithium
Charger	: 2A
Motors	: 24V 150W x 2
Max Gradient	: 6 degrees
Wheels	: Front 7" Solid, Rear 12" Solid