

MOBILITY SCOOTERS



MOVING LIFE ATTO FOLDABLE MOBILITY SCOOTER (CLASSIC)

- LTA-compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Stylish contemporary design
- Easy to fold into suitcase size and easy to Unfold
- Adjustable seat with 3 different heights
- Easily adjust height and angle of tiller
- Lightweight aluminium chassis
- Flat-free tyres for worry-free travel.

SPECIFICATIONS

Size Folded	: 38.9 x 41.9 x 71(cm) / 15.3" x 16.5" x 28"
Size Unfolded	: 120 x 56 x 71 (cm) / 47" x 22" x 35"
Weight Front	: 12 kg / 26.4 lbs
Weight Rear	: 17 kg / 37.5 lbs
Assembled Weight	: 29 kg / 64 lbs (w/o battery)
Battery Weight	: 1.85 kg / 4.1 lbs
Total Weight	: 30.85 kg / 68 lbs
Speed	: 6.4 km/h / 4 mph
Capacity	: Maximum user weight 265 lbs / 120 kg
Range	: Up to 20 km/12.5 miles (continuous driving)
Battery	: 48 V, 5.2 AH lithium-ion
Battery Size	: 249.6 wh
Charger	: International 110-240V AC charger. Charge on-board and off-board.
Cruise Control (Throttle Not Needed)	: No
Automatic Slowdown	: Slopes
Front Disk Brake for Extra Safety	: No
Display	: LCD / Programmable
Tires	: Airless PU, Flat free tires
Antitheft Wireless Locking FOB Key	: No
Embedded Front Headlights	: No
Seat and Backrest / grips	: PU Seat / Grips
Rear Light Lifting Handle	: No
Eco Mode	: No
Left Handed Operation	: Optional