

# MOBILITY SCOOTERS

## EFOLDI LIGHTWEIGHT FOLDING MOBILITY SCOOTER

- LTA-compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Weighs as light as 12 kg with rear wheels and battery removed.
- Can be pulled along like a trolley
- When fully folded, it is lighter and smaller than a mid-sized luggage bag.
- Travel range of up to 16 km



### SPECIFICATIONS

Folded size (L/W/H)	: 65x55x37(cm),26x22x15 (in)
Gross Weight with Battery	: 117x55x94(cm),46x22x37(in)
Seat height	: 51cm (20inch)
Footrest height	: 12cm (4.7inch)
Front wheel	: 20cm (8inch) Solid
Rear wheel	: 22.5cm (9inch) Solid
Weight with battery	: 17kg (38lbs)
Weight without battery	: 15kg (33lbs)
Weight without battery & rear wheels	: 12kg (26lbs)
Maximum Speed	: 6.4 km/h (4 mph)
Maxcarry capacity	: 120kg (265lbs)
Maximum Incline	: 1:10 Gradient (6 degrees)
Turning Radius	: 100 cm (40inch)
Maximum kerb height	: 4cm (1.5inch)
Motor	: 180W, 24VDC Hub motor
Battery	: 24VDC,10Ah Lithium Ion, 240Wh
Charger	: 110-220VAC, 2A
Maximum range	: 10miles (16km)
Brake	: Electromagnetic