

MOBILITY SCOOTERS

PHOENIX HD 3-WHEELED PORTABLE MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- 3 Wheel configuration allows scooter to nimbly negotiate indoor obstacles
- Comfortable, height adjustable seat with flip back, adjustable width arms
- Flat-free tyres for worry-free travel.
- Adjustable angle tiller
- Comes with easy to remove, quick connect battery pack that can be charged separately from the scooter.



SPECIFICATIONS

Dimensions (l x w)	: 38.5" x 21.7" / 98 cm x 55 cm
Base weight	: 66.6 lbs / 30.3 kg
Battery weight	: 30.6 lbs / 13.8 kg
Seat weight	: 23.8 lbs / 10.8 kg
Gross weight	: 121lbs/55kg
Front Wheel Dimension	: 8"x2"/20cm x5cm SOLID
Rear Wheel Dimension	: 9"x3"/22cm x8cm SOLID
Battery (included)	: 20 AH, 12V X 2
Charger	: 2 A Off Board
Maximum speed	: 6 kmh
Maximum capacity	: 350 lbs / 158 kg
Cruising range	: Upto12.5miles/20km
Turning radius	: 43"/110cm
Seat Dimension	: 18"/46cm (Width), 17.5" / 44 cm (Depth)
Backrest height	: 15" / 38 cm
Seat-to-deck height	: 13.5"-15.5"/34-39 cm Adjustable
Seat-to-floor height	: 20"-22"/51 - 56 cm Adjustble
Armrest Dimensions	: 9" / 23 cm (Height), 12.5" / 32 cm (Length)
Armrest-to-armrest width	: 19"-23"/48-54 cm Adjustable
Climbing angle	: 8°
Controller	: PG 70A