

MOBILITY SCOOTERS



SOLAX GENIE AUTOMATIC FOLDING MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Scooter can be easily folded within seconds using a remote control.
- 15 km maximum range with standard 24V 10AH Lithium battery. Range can be increased to 30 km with optional spare battery.
- Flat-free tyres for worry-free travel.
- Wrap-around delta handles make driving easy, even for those with who can only use one hand.
- Adjustable angle tiller allows the user to position the controller into the most ideal driving position.

SPECIFICATIONS

Unfolded Dimensions (mm)	: 930 (L) X 450 (W) X 857 (H)
Folded Dimensions (mm)	: 365 (L) X 450 (W) X 580 (H)
Seat Height	: 560 mm
Tire (Front wheel)	: 6 inches(solid wheel)
Tire (Rear wheel)	: 7 inches (solid wheel)
Maximum Speed	: 6 kmh/3.75 mph
Safe Gradient	: Up to 15 Degrees
Maximum Gradient	: 8 Degrees
Range	: 15 km/9.32 miles
Turning Circle	: 1200 mm
Total Weight	: 25 kg (55lbs) (Battery included)
Brake	: Intelligent regenerative electromagnetic brake
Drive System	: Rear wheel drive
Handlebar	: 3 positions
Maximum Capacity	: 120 kg
Ground Clearance	: 115 mm (4.5 inches)
Battery	: 24 V/10 AH (Lithium Battery)
Motor	: 120 W
Battery Charger	: DC 24V/2 A off-board (automatic protection charger)