

MOBILITY SCOOTERS

MOBIFREE FOLDING MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Scooter can be easily folded within seconds.
- 15 km maximum range
- Flat-free tyres for worry-free travel.
- Wrap-around delta handles make driving easy, even for those with who can only use one hand.
- Adjustable angle tiller allows the user to position the controller into the most ideal driving position.



SPECIFICATIONS

Size (Unfolded) (mm)	: 969 (L) X 530 (W) X 950 (H)
Folded Dimensions (mm)	: 343 (L) X 530 (W) X 950 (H)
Front Wheels	: 6" x 2" Solid Tyres
Rear Wheels	: 7" x 2" Solid Tyres
Seat Size (cm)	: 42 (W) x 40 (D) x 32 (H)
Seat-to-Floor Height	: 53 cm
Max Speed	: 6 km/h
Max Gradient	: 12 Degree
Max Range	: Up to 15 km
Turning Radius	: 1.4 m
Weight excluding battery	: 19.2 kg
Weight including battery	: 21.0 kg
Maximum Capacity	: 125 kg
Battery	: 24V, 10H (Lithium Battery)
Charger	: DC 24V / 2A with Safety Mark