

MOBILITY SCOOTERS



MOJO AUTOMATIC FOLDING MOBILITY SCOOTER

- LTA compliant personal mobility aid (PMA).
- Charger comes with Singapore Safety Mark.
- Scooter can be easily folded within seconds using a remote control.
- 15 km maximum range provided by 2 pieces of 24V, 5Ah lithium batteries (7.5 km each)
- Flat-free tyres for worry-free travel.
- Wrap-around delta handles make driving easy, even for those with who can only use one hand.
- Adjustable angle tiller allows the user to position the controller into the most ideal driving position.

SPECIFICATIONS

Unfolded Dimensions (mm)	: 1010 (L) X 450 (W) X 935 (H)
Folded Dimensions (mm)	: 490 (L) X 450 (W) X 690 (H)
Tire	: 7 inch, solid honeycomb
Maximum Speed	: 6 km/h / 3.75 mph
Maximum Gradient	: 12°
Range	: 15 km / 9.32 miles
Turning Radius	: 1400 mm
Weight - Front Section	: 19 kg / 41.9 lbs
Weight - Rear Section	: 9 kg / 19.8 lbs
Weight - Batteries	: 2.5 kg / 5.5 lbs
Total Weight	: 30.5 kg / 67.2 lbs (With Battery)
Brake	: Intelligent regenerative electromagnetic brake
Drive System	: Rear wheel drive
Maximum Capacity	: 125 kg / 275 lbs
Battery	: 24V, 5AH(Lithium Battery)x2 pcs
Motor	: 120 W
Battery Charger	: DC 24V / 2A