

Meet the iGo **EZ Chair** Family



EZ Chair Lux

Top of the range electric wheelchair, that is powerful enough to go anywhere, fold up into a boot of a car and travel on an airline. Carries heavier weights and has provision with a wider seat, super powerful to drive on any terrain outdoors.

EZ Chair Base

Value folding scooter that fold in seconds and fits into most boots. Use indoors anywhere and level outdoor terrain.

EZ Chair Mini

Quick folding scooter, built for shopping and mall visits. Lightweight and its small size makes this scooter unique.

SPECIFICATIONS

	EZ Chair Lux	EZ Chair Base	EZ Chair Mini
Motor type	Brushless	Brushed	Brushless
Motor power	270	250	250
Battery size	6AH x2 =12AH	12AH	10AH
scooter size folded	75x60x34	75x63x37	73x71x26
Scooter size built	100x94x59	110x96x63	100x95x60
Scooter weight	28	29	20
Seat width	52	45	45
Max Weight can carry	120	90	80
Airline approved*	Yes	No	No
Climbing angle	20 DEG	13 DEG	10 DEG
Charge time	6HRS	6HRS	4HRS
Waranty: Battery	6 months	6 months	6 months
Waranty: Electronics	2 year	1 year	1 year
Warranty: working parts	2 year	1 year	1 year