

# iGo Mobility

SUPERIOR DESIGNS | FIRST CLASS LUXURY



COMPACT  
PORTABLE



LARGE  
PORTABLE



FULL-SIZE  
RUGGED LUXURY

# iGo Mobility



# iGo Compact

Lite GB106



XL GB116

COMPACT PORTABLE

GB146

The iGo Compact Series is perfect for those individuals seeking a fast and convenient way to travel. Do you have a long day of walking planned? Not a problem for the iGo Compact, simply assemble/disassemble and get ready for a day full of cruising light.

## iGo Compact

- 300 lb. weight capacity
- Delta Tiller standard
- Disassembles into 5 lightweight pieces
- Battery charging port conveniently located below the control panel
- Top speed 4.5 mph
- Wireless and easy disassembly!
- Comfortable 17"x16" seat standard
- Two sets of interchangeable shroud panels **INCLUDED FREE:** Red and Blue
- Under carpet troubleshooting guide included

## iGo Compact Lite

- 250 lb. weight capacity
- Delta Tiller standard
- Disassembles into 5 lightweight pieces
- Top speed 4 mph
- Three sets of interchangeable shroud panels **INCLUDED FREE:** Red, Blue and Sahara

# COMPANION



FULL-SIZE LUXURY

The Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- LiquiCell™ comfort equipped seat
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Backlit battery gauge
- Black, non-marking, low profile tires on sturdy steel rims
- Comfortable, fully adjustable padded armrests
- Under-carpet troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC340

Companion GC440

Companion GC240

# iGo Mobile



GL110

LARGE PORTABLE

GL140

The unique design of the iGo Mobile scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble when needed, and maneuver smoothly indoors and out. Both iGo Mobile models feature:

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Adjustable tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat
- Adjustable armrests
- 18Ah batteries
- Completely wireless disassembly
- 9" front and rear tires
- Low-profile, black tires on high-strength steel rims with stylish wheel covers





## FULL-SIZE RUGGED LUXURY

When you want a rugged, heavy duty scooters to tackle the challenging terrain of the great outdoors, look no further than the **iGo Outdoor** and **Avenger™**. Both feature oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.

### **Avenger**

The Avenger model features:

- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- 18-mile operating range

## iGo Outrider

The **iGo Outrider** model features:

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge
- Black, low profile, non-marketing solid tires with chrome hubcaps



## AVAILABLE ACCESSORIES

\*Only available on Companion models.



Comfort Ride\*



Crutch Holder



Rear Basket



Oxygen Holder



Walker Holder



Safety Flag



Cup Holder



Pack'n'Go



Power Elevating Seat\*



Quad Cane Holder

# SPECIFICATIONS\*

<i>Scooter Model</i>	<i>Length</i>	<i>Width</i>	<i>Tires</i>	<i>Front Tire</i>	<i>Rear Tire</i>	<i>Weight Capacity (lbs.)</i>	<i>Top Speed (mph)</i>	<i>Ground Clearance</i>	<i>Operating Range** (- miles)</i>	<i>Turning Radius</i>	<i>Batteries (2)</i>	<i>Weight w/ Batteries</i>	<i>Front Half Weight</i>	<i>Rear Half Weight</i>	<i>Seat Weight w/ Arms</i>	<i>Heaviest Piece</i>	<i>Freewheel Mode</i>	<i>Electronic Control</i>	<i>Charger</i>
GB-106	36.5"	21.5"	Solid Low Profile	8"	8"	250	4	1.75"	8	31"	2-12Ah	96	29	26	22	29	Yes	Yes	Off-board
GB-116	40.5"	21.5"	Solid Low Profile	8"	8"	300	4.5	2"	7/10 <sup>1</sup>	37"	2-12 Ah/ 2-18 Ah	102	33	25	22	33	Yes	Yes	Off-board
GB-146	42"	21.5"	Solid Low Profile	8"	8"	300	4.5	2"	7.5/11 <sup>1</sup>	48.5"	2-12 Ah/ 2-18 Ah	109	41	25	22	41	Yes	Yes	Off-board
GL-110	41.5"	22"	Solid Low Profile	9"	9"	300	5	3.75"	9.5	39"	2-18Ah	119	44	28	22	44	Yes	Yes	Off-board
GL-140	43.25"	22"	Solid Low Profile	9"	9"	300	5	3.75"	9.5	48"	2-18Ah	127	52	28	22	52	Yes	Yes	Off-board
GC-240	42.5"	23.5"	Solid Low Profile	9"	9"	300	4.5	3.75"	15.5	39"	2-U1	166	46	48	28	48	Yes	Yes	Off-board
GC-340	45.5"	23.5"	Solid Low Profile	9"	10"	350	4.5	4.5"	16.5	40.5"	2-U1	168	48	48	28	48	Yes	Yes	Off-board
GC-440	47.5"	23.5"	Solid Low Profile	10"	10"	350	4.5	4.5"	15	53.5"	2-U1	183	63	48	28	63	Yes	Yes	Off-board
GA-541	57.5"	24.5"	Pneumatic	13"	13"	500	7	4"	18	67"	2-NF22	265	68	68	55	68	Yes	Yes	Off-board
GR-575	57"	28"	Solid Low Profile	13"	13"	400	8	3.5"	18	76.5"	2-75AH	314	Base 158	N/A	50	160	Yes	Yes	Off-board

\* ALL WEIGHT SHOWN IN LBS.

\*\*Battery range will vary due to rider weight, drive surface and drive slope. Range was tested at maximum weight.

<sup>1</sup>. Range tested at 200 lbs.

**Brakes:** All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems. **Rear wheel drive:** All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.