

Meet the iGo **EZ Chair** Family



EZ Chair Lux

Top of the range electric wheelchair, that is powerful enough to go anywhere, fold up into a boot of a car and travel on an airline.

Carries heavier weights and has provision with a wider seat, super powerful to drive on any terrain outdoors.

EZ Chair Base

Value folding scooter that folds in seconds and fits into most boots. Use indoors anywhere and level outdoor terrain.

EZ Chair Light

Quick folding scooter, built for shopping and mall visits. Lightweight and its small size makes makes this scooter unique.

SPECIFICATIONS

	EZ Chair Lux	EZ Chair Base	EZ Chair Light
Motor type	Brushless	Brushed	Brushless
Motor power	270	250	180
Battery size	6AH x2 =12AH	12AH	12AH
scooter size folded	75x60x34	75x63x37	73x71x26
Scooter size built	91x61x100	95x57x100	58x33x70
Scooter weight	28	29	15
Seat width	52	45	40
Max Weight can carry	120	90	80
Airline approved*	Yes	No	No
Climbing angle	20 DEG	13 DEG	10 DEG
Charge time	6HRS	6HRS	6HRS
Waranty: Battery	6 months	6 months	6 months
Waranty: Electronics	2 year	1 year	1 year
Warranty: working parts	2 year	1 year	1 year