

ScopeValet™

Endoscopy Care Products

ScopeValet™

Disposable Biopsy Valves

Item No. 345BVO - Olympus and Fujinon

Item No. 345BVP - Pentax

ScopeValet™ DISPOSABLE BIOPSY VALVES are non-sterile and used to cover the openings to the biopsy/suction channel inlets of Olympus®, Fujinon® or Pentax® gastrointestinal endoscopes.

The biopsy valves provide access for endoscopic device passage and exchange, help to maintain insufflation, and minimize the leakage of bioburden from the biopsy port throughout the gastrointestinal endoscopic procedure. The disposable biopsy valves are designed for one time use which eliminates the need for manual cleaning and reprocessing. Most important the product reduces cross contamination helping to prevent potential patient/staff safety risks.



Olympus® and Fujinon®
Disposable
Biopsy Valves



Pentax®
Disposable
Biopsy Valves

Features

- Minimal "spitting" or leakage for general GI procedures
- Maintains scope pressure while passing devices
- Helps reduce cross-contamination
- Eliminates the need for reprocessing dirty valves
- Color-coded and designed to work with the major endoscopes in the market:
RED to fit Olympus® & Fujinon® scopes; WHITE to fit Pentax® scopes.
- Latex free
- Disposable, one time use

Contraindication:

This device is for one time use only. The reuse of single use devices might present risks to the safety of individuals, due to possible uncontrolled contamination.

*Olympus® is a registered trademark of Olympus Optical Co., Ltd.

*Fujinon® is a registered trademark of Fujinon Inc.

*Pentax® is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha.