

Ruhof InstruFlush™

Item No. 345IF

Independent Flushing System for Effective Channel/Lumen Pre-Cleaning & Flushing

Ruhof InstruFlush™ is an effective instrument and scope channel/lumen pre-cleaning and flushing device which increases the productivity of syringe flushing by 74.9%. The Ruhof InstruFlush™ provides a pulse feature for hard-to-dislodge soils and works very well with low-foaming enzymatic chemistries and detergents. The pump can be used with or without the Ruhof InstruStation™.

Features

- Consistently meets ANSI/AAMI ST79: 2017 recommendations for cleaning, flushing, and rinsing
- Offers components (measurements on reservoirs and in sink; sink thermometers) that facilitate compliance to cleaning solution and medical device manufacturers' IFU
- Flushes water or cleaning solution
- Can be used with or without the Ruhof InstruStation™
- Up to three pumps can be mounted at one time for optimized capacity and productivity (Three pumps can flush all lumens in a complex or robotic tray all at once, in one minute)
- Can dispense up to 1.2 liters of cleaning solution in ONE MINUTE
- Has an Impact Cycle to help remove soil from both ends
- Automatically flushes up to 1,000ml of fluid in one minute; three speed controls enable you to flush a variety of instruments; assures flushing consistency
- Three lumens or channels can be flushed at one time by one pump
- Integrated pressure relief protects scopes
- Includes accessories to wall mount one, two or three pumps in a row
- Greatly reduces risk of repetitive motion injuries
- Frees up technician to perform other tasks



Pump can be wall mounted, one, two or three pumps in a row, to promote a clean work area.



InstruFlush™ will pre-clean both flexible and rigid scopes, large and small.



InstruFlush™ can flush robotic instruments and can be used with or without Ruhof InstruStation™.

PLEASE REFER TO InstruFlush IFU for complete instructions on how to use Ruhof's InstruFlush Independent Flushing System.