

New from ADC

Introducing Our Most Eye-catching Stethoscope Ever:


The Adscope® 603 Iridescent

Our most popular clinician scope, the Adscope 603, is now available with a rainbow-inspired, iridescent chestpiece. This luminous finish can be paired with our Metallic Caribbean or Metallic Raspberry tubing, and is fitted with a clinician headset in a coordinating golden tone. A tactical version is also available.

The Adscope 603 offers:

- Combination diaphragm/bell adult chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- Ultra-sensitive diaphragm (1 3/4" diameter) with AFD technology and an extra-large bell (1 3/8" diameter)
- Clinician headset boasts reinforcing yoke molded into the flexible 22" PVC tubing and aluminum binaurals fixed at 15° angle for greater comfort
- Includes accessory kit with two extra pairs of spare Adsoft™ Plus eartips and scope ID tag
- Weighs 5.8 oz. (164 g), measures 31" long
- Assembled, inspected, and packaged in the USA
- Lifetime warranty and free parts forever

Visit www.adctoday.com/products/603
for additional details.



603MCA
Iridescent Metallic Caribbean
MSRP: \$67.99



603MRS
Iridescent Metallic Raspberry
MSRP: \$67.99



603IST
Iridescent Tactical
MSRP: \$67.99

55 Commerce Drive
Hauppauge, New York 11788
Tel: 1-800-ADC-2670
Fax: 631-273-9659
www.adctoday.com


AMERICAN DIAGNOSTIC CORPORATION