

UNLEASH THE POWER OF SINGLE-USE

ABOUT AMBU

Since 1937, our desire to make life easier for doctors and create innovative products that save lives has driven Ambu. We continue to build on our heritage by driving breakthrough ideas forward and turning them into efficient healthcare solutions. It has brought us to where we are today – the world's most innovative single-use endoscopy company.



Ambu® Product Portfolio

Ambu

ANESTHESIOLOGY AND RESPIRATORY

Our latest innovations for anesthesiologists will help save lives and improve patient care. This full portfolio of quality anesthesia products crosses over a variety of applications, from everyday usage to exceptional situations.




Masks

The King Mask is Ambu's premium single-use disposable anesthetic face mask. The anatomically shaped mask features a thin cushion that is tacky to the touch which supports a tight seal, the top of the cushion is ribbed for a no-slip grip for the clinician.



Breathing Circuits

The innovative design of the King F family of circuits warm and humidify the gas delivered to the patient, and can reduce the number of SKUs a hospital must stock. The King Flex2, F Flex, and F Flex2 are the only single-limb circuits to utilize flex tubing, which can be expanded to multiple lengths, further reducing the hospital's SKU count.



Pain Pumps

The advanced regulation technology of the ACTion disposable infusion pump provides consistent pain relief for the patient and allows anesthesia providers to address multiple clinical needs with one pump, reducing cost and inventory.



Breathing Bags

Ambu's King Breathing Bags are designed with the clinician in mind, providing an excellent feel for lung compliance during ventilation. They come in standard green but also in a range of colors for pediatric patients. They are thin, but durable and latex free. The hourglass-shaped breathing bag efficiently delivers required gas volume.



Ambu® SPUR® II Resuscitators

Ever since the invention of the world's first resuscitator in 1956 – the now iconic Ambu bag – Ambu has become the epitome of reliability and convenience.



Laryngeal Tubes

The Ambu® King LTS-D™ Laryngeal Tube is a single-use alternative airway device that is simple to use and provides superior patient ventilation. Designed with a straightened, beveled distal tip that helps in directing the airway posterior to the larynx and into the upper esophagus, it enables passage of the gastric tube through a separate channel. Its unique configuration means minimal risk of the device entering a patient's trachea. The Ambu device is a powerful tool in improving out-of-hospital cardiac arrest outcome. The King LTS-D Laryngeal Tubes are available in seven sizes ranging from 0 to 5, and in kits which include syringe and lube.

Laryngeal Masks

The AuraGain is Ambu's 2nd-generation laryngeal mask. It satisfies three fundamental airway management needs by integrating gastric access and intubation capability in an anatomically curved single-use device that facilitates rapid establishment of a patent airway.


**COMBINED
110 YEARS
OF INNOVATION
IN ANESTHESIA**





PATIENT MONITORING AND DIAGNOSTICS

Our patient monitoring and diagnostics product line includes innovative solutions for your cardiology, neurology, and intraoperative monitoring needs.




ASK US ABOUT OUR SINGLE PATIENT PACKAGING OPTIONS FOR CUP ELECTRODES

Neurology


Ambu has worked for decades with neurophysiology experts to co-develop tools that meet the quality you demand.



Cup Electrodes



EMG Guided Injections




Surface Electrodes



EMG Needles

Intraoperative Monitoring


Protect patients during surgery by monitoring their nervous system with our line of surface electrodes, probes, and subdermal needles.



Corkscrew Electrodes



Surface Electrodes



Subdermal Needles



Probes

Cardiology

Since the first single-use Ambu® BlueSensor™ ECG electrodes were introduced in 1976, we have delivered quality solutions matched to your specific needs. Our full portfolio, with its wide range of sizes, shapes, gel types and backing materials, provides devices for every situation.



Neonatal



Pediatric




Adult



AIRWAY MANAGEMENT

The inventor of the world's first self-inflating resuscitation bag offers a full portfolio of products for airway management professionals to help save lives and improve patient care. Explore the clinical evidence and learn more about our products and resources.




aCart™ Compact

Ambu's new aCart™ Compact is a modern, mobile workstation for healthcare workers who perform endoscopic procedures. It provides easy storage for our single-use product portfolio in a small footprint. The aCart is convenient, portable, versatile, robust, and hygienic.



Video Laryngoscopes

King Vision is designed for indirect visualization of the vocal cords which minimize the need for tissue manipulation and improves the ability to see anterior airways. Such reliable performance is a necessity, not an option. King Vision comes with the choice of two blade types that allow users to keep their preferred technique in terms of tube delivery, offering a unique flexibility.



Fuji UNIBLOCKER™

UNIBLOCKER is the only bronchial blocker that incorporates a steel mesh polyurethane covered shaft for direct control of the blocker into the targeted bronchus. This control allows for simple and intuitive placement of the blocker.

SINGLE-USE ENDOSCOPY

Ambu created the world's first single-use flexible endoscope in 2009. We now leverage more than a decade of single-use flexible endoscopy experience in multiple clinical areas, with more to come in the near future. Join the world leader in single-use endoscopy on this journey.




Pulmonology

Sterile straight from the pack, Ambu bronchoscopes are part of an integrated solution that helps safety, workflow and productivity. With aView™ 2 Advance HD monitor, capture images and video. We're extending our solution to the bronchoscopy suite with our 5th-generation scope. Be equipped as regulatory agencies increasingly recommend single-use.



aScope™ 5 Broncho



aScope™ 4 Broncho



aScope™ 4 Broncho Sampler Set



aView 2 Advance

Gastroenterology

Single-use gastrointestinal endoscopes remove the risk of patient cross-contamination. Ambu believes that sterile single-use gastrointestinal endoscopy solutions should perform like conventional ones, only with the added advantages of increased availability and improved patient safety.



aScope™ Duodeno



aScope™ Gastro




aBox™ 2

Urology

Ambu provides a single-use solution for urology to increase your productivity and simplify your workflow. The aScope™ 4 Cysto is a sterile single-use, flexible cystoscope containing an all-in-one HD display with excellent imaging, smooth maneuverability, and patient-ready out of the package.



aScope™ 4 Cysto



aView 2 Advance

ENT

Ambu offers single-use rhinolaryngoscopes with high-quality imaging, multiple connectivity options for outpatient clinics and easy portability for inpatient consults. Convenient for fiberoptic endoscopic evaluation of swallowing (FEES) procedures, our scopes offer frame-by-frame video replay and, soon, audio recording.



aScope™ 4 RhinoLaryngo



aView 2 Advance



Ambu FOREVER
FORWARD

Questions? Contact your Ambu representative at
800-262-8462 or visit **www.ambuUSA.com**

Ambu Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

US: RX only