



Easy TBI
for the
Diabetic
Foot

dopplex®

Ankle and Toe Pressure Kit

Simple and Rapid Assessment of the Diabetic Foot

A comprehensive kit for vascular assessment of the diabetic foot and lymphatic limb.

The ankle and toe pressure (ATP) kit is a unique combination of the DMX Digital Doppler, probes and accessories providing the vascular, diabetic or lymphatic Specialist with a system to aid the assessment of arterial disease and neuropathy.

Features:

- Complete kit for easy measurement of the Toe Brachial Index (TBI) on the diabetic and lymphatic foot
- DMX Doppler features a high-resolution color display of the vascular waveform
- APPG probe for easy measurement of TBI
- DMX Doppler features a Digital Dynamic Noise Reduction system to reduce background noise
- Ability to connect to optional Dopplex Reporter 5 reporting software to store and view results electronically

Included in the box:

- DMX Digital Doppler
- Medical Grade Recharging Kit
- EZ8XS Wide beam 8 MHz Doppler Probe
- Toe cuffs and inflator
- APPG probe for toe
- Adhesive strips for APPG probe (pack of 240)
- Arm/Ankle Cuffs (standard and large; latex free)
- Trigger operated Sphyg
- Neuropen and box of tips
- Educational CD
- Large Carry bag

Our arterial plethysmography (APPG) probe easily attaches to the DMX Doppler Unit by simply swapping your Doppler probe with its adaptor.

Configurations

ATP-USA Ankle Toe Pressure Kit

Warranty

3 years on Doppler Unit

1 year on probe, adaptors, cable and accessories

Functionality

Complete kit for assessment of Arterial Disease or Neuropathy of the diabetic and lymphatic foot

Measuring Toe Brachial Index (TBI)

Electronically store and print results via the optional Dopplex Reporter 5 (DR5) software



Contact:

Arjo, Inc. 2349 W. Lake St., Suite 250, Addison, IL 60101
Phone: +1 (800) 323 1245 Fax: +1 888 594 2756
Web: www.huntleigh-healthcare.us Email: Diagnosics.USADD@arjo.com

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY