



- IV Start Kits
- Sterile Peripheral IV Catheter Insertion Kits (Ster-Assist™)
- Extension Set/ IV Securement and Dressing Kits
- Dressing Change Kits

B. Braun IV Accessory Kits

Convenient and Compliant IV Care
at Your Fingertips

Find out more about our standard and custom kits.

Call your B. Braun Representative at 1-800-227-2862

1-800-227-2862 | www.BBraunUSA.com

Rx Only. ©2016 B. Braun Medical Inc., Bethlehem, PA. All Rights Reserved.
ChloroPrep is a registered trademark of Cardinal Health, Inc., 3M and Tegaderm are a registered trademarks of 3M,
CIV-Flex is a trademark of CIVCO Medical Solutions.

VAD82_1/16_BB

B. Braun provides all the necessary components you need in convenient kits.

Kit components are designed to help reduce catheter-related complications, and maintain compliance when placing a peripheral IV line or changing a central line dressing.

IV Start Kits

Quick and easy access to components needed to complete a peripheral IV insertion.



Kits include: tourniquet, gauze, sterile transparent dressing, and tape and are available with other popular items, such as:

- B. Braun's needleless extension sets with Caresite® Luer Access Device
- ChloroPrep®
- 3M® Tegaderm®

Customize your kit by mixing and matching available components.

Sterile Peripheral IV Catheter Insertion Kits (Ster-Assist™)

Comprehensive kits that support sterile and ultrasound guided procedures.



Kits include all the necessary components required to support an ultrasound guided PIV insertion, and sterile technique.

- B. Braun's longer length Introcan Safety® IV catheters
- B. Braun's Caresite needle-free extension set
- CIV-Flex™ transducer cover and ultrasound gel
- 3M Tegaderm IV advanced securement dressing 1683

No customizable options available for these kits.

Extension Set/IV Securement and Dressing Kits

Needleless IV access, securement, and protection for peripheral IV catheter site care.



Kits include the key components needed to complete the peripheral IV insertion process.

- B. Braun's Caresite needle-free extension sets
- IV securement dressing
- Alcohol skin prep

Complies with CDC¹ and INS² standards and guidelines for peripheral IV catheter.

Customize your kit by mixing and matching available components.

1. CDC Guidelines for the Prevention of Intravascular Catheter-related Infections, 2011.
2. INS Standards of Practice 2011, Std. 36, Vascular Access Device Stabilization

Dressing Change Kits

Sterile components designed to safeguard your patient from infection during site dressing changes.



Kits include:

- CSR wrap or drape
- Mask and gloves
- ChloroPrep® or alcohol swabsticks
- Gauze and tape
- Sterile transparent dressing

Customize your kit by mixing and matching available components.