


WE COULD MAKE THEM THIS WAY... BUT WE DON'T!

Bovie®

BOVIE MEDICAL CORPORATION • 5115 Ulmerton Road, Clearwater, FL 33760 • USA • Ph 1-800-537-2790 • Sales@BovieMed.com • www.BovieMedical.com



We could make them this way...but we don't! We make them better!

1. We could use heat shrink instead of a double molding process...But we don't.
2. We could leave off the molded ribs for easy insertion and removal from the pencil...But we don't.
3. We could leave out the patented protective tray for the Tungsten loop...But we don't.
4. We could make them in China...But we don't. We make them in the U.S.A.
5. We could skip the 100% inspection, which includes running 4500 volts through the electrode for 30 seconds...But we don't.
6. We could charge a lot more for our loop electrodes...But we don't.



There are many steps we could leave out of our process for making a great loop electrode...but we don't!

LIST OF ELECTRODE STYLES AND CONVERSION CHART...

BOVIE / AARON ELECTRODES			COOPER	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CIRCON	UTAH	MILTEX
Ball Electrodes	ES06	Extended 3mm Ball, 5/box	B0312		909038	60-5181-102	E1563	6156901	DBL311	31-5000
	ES07	Extended 5mm Ball, 5/box	B0512		909003	60-5181-103	E1564	6157901	DBL511	31-5002
LLETZ Loop Electrodes	ES08	5mm x 5mm Loop, 5/box						6133901	DLP-T11	
	ES09	10mm x 10mm Loop, 5/box	R1010		909007		E1559	6134901	DLP-S11	31-5056
	ES10	15mm x 8mm Loop, 5/box			909011			6135901		31-5072
	ES11	15mm x 10mm Loop, 5/box				60-5180-102		6136901		31-5074
	ES12	20mm x 8mm Loop, 5/box	R2008		909013			6137901		31-5090
	ES13	20mm x 15mm Loop, 5/box	R2015		909009			6137901		31-5096
	ES31	20mm x 20mm Loop, 5/box	R2020							31-5098
	ES42	20mm x 12mm Loop, 5/box					E1561			31-5094
	ES43	15mm x 15mm Loop, 5/box	R1515							
	ES44	15mm x 5mm Loop, 5/box	R1505		909134					
	ES45	13mm x 13mm Loop, 5/box								
	ES46	10mm x 5mm Loop, 5/box								
	ES47	25mm x 10mm Loop, 5/box	R2510							
	ES49	20mm x 10mm Loop, 5/box	R2010		909132	60-5180-104				
	ES51	13mm x 8mm Loop, 5/box				60-5180-101*				
	ES52	10mm x 8mm Loop, 5/box								31-5054
ES53	15mm x 12mm Loop, 5/box					E1560				
LLETZ Square Electrodes	ES14	5mm x 5mm Square, 5/box					E1565	6153901		
	ES15	10mm x 4mm Square, 5/box	S1004				E1566	6154901	DLP-SQ2	
	ES16	10mm x 8mm Square, 5/box			909017		E1567	6155901	DLP-SQ1	
	ES41	5mm x 10mm Square, 5/box								

* shows 1mm difference in size

ELECTRODE SHAPES AND SIZES...

