



ES22-8 ES23-8 ES24-8 ES25-8 ES26-8

DISPOSABLE SHORT TUNGSTEN .008 WIRE LOOPS

- ES22-8 — 1/4" short shaft loop, sterile, .008" wire - 5/box
- ES23-8 — 5/16" short shaft loop, sterile, .008" wire - 5/box
- ES24-8 — 3/8" short shaft loop, sterile, .008" wire - 5/box
- ES25-8 — 1/2" short shaft loop, sterile, .008" wire - 5/box
- ES26-8 — 3/4" short shaft loop, sterile, .008" wire - 5/box



ES22-R ES23-R ES24-R ES25-R ES26-R

REUSABLE SHORT TUNGSTEN LOOPS


- ES22R — 1/4" short shaft loop, non-sterile - 1/box
- ES23R — 5/16" short shaft loop, non-sterile - 1/box
- ES24R — 3/8" short shaft loop, non-sterile - 1/box
- ES25R — 1/2" short shaft loop, non-sterile - 1/box
- ES26R — 3/4" short shaft loop, non-sterile - 1/box

...FROM A NAME YOU CAN TRUST

All Aaron® blade, ball, needle, loop, and square disposable electrodes use the standard 3/32" diameter stainless steel shafts. These electrodes, like all Aaron® electrodes, feature safety grip insulators, combining patient and user safety with easy insertion into and removal from the surgical pencil.



Many Electrode Sizes Available



The loop and square electrodes feature tungsten wire for superior shape and wire integrity throughout the excision procedure, while minimizing thermal damage and yielding excellent histopathology specimens. Each molded loop and square electrode includes our patented tray to protect the loops from bending during shipping and handling.



Ball electrodes in 3mm and 5mm sizes are available for fulguration and desiccation, while needle electrodes are available for precise, pinpoint cutting and coagulation.

Distributed By:



7100 30th Aveue North • St. Petersburg, FL 33710-2902
 U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640
 Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144
 www.boviemed.com • sales@boviemed.com

List of Electrode Styles and Conversion Chart...




			STRYKER/ COLORADO	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CIRCON	UTAH	MILTEX
SuperCut™ Tungsten Needles	ES60	3CM Super Fine Needle, 5/box	N101A							
	ES61	3CM Modified Super Fine Needle, 5/box	N103A				E1651			
	ES62	4.5CM Super Fine Needle, 5/box	N102A							
	ES63	4.5CM Modified Super Fine Needle, 5/box	N104A							
			COOPER	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CIRCON	UTAH	MILTEX
Needle Electrodes	ES02	Standard Needle, 25/box			909135	275-1502/9113	E1552			
	ES03	Extended Needle, 25/box				275-1504/9116	E1552-6			
	ES38	Modified Needle, 25/box								
	ES40	Extended Modified Needle, 25/box	N2512	9006240	909005	138025		6162901		
Dermal Loops	ES22-8	1/4" Short Shaft Loop, 5/box		9006406		138011				
	ES23-8	5/16" Short Shaft Loop, 5/box								
	ES24-8	3/8" Short Shaft Loop, 5/box		9006407		138010				
	ES25-8	1/2" Short Shaft Loop, 5/box		9006400						
	ES26-8	3/4" Short Shaft Loop, 5/box				138012				
Dermal Tips	A804	Sharp Dermal Tip, Non-Sterile, 100/box								
	A805	Sharp Dermal Tip, Sterile, 50/box								
	A806	Blunt Dermal Tip, Non-Sterile, 100/box								
	A807	Blunt Dermal Tip, Sterile, 50/box								
Blade Electrodes	ES01	Standard Blade, 50/box				275-1500/9111	E1551G			
	ES04	Extended Blade, 25/box				275-1503/9114	E1551-6			
	ES37	Modified Blade, 25/box				275-1524	E1553			
	ES39	Extended Modified Blade 25/box				275-1525				
	ES18	Angled Blade, 5/box								
Ball Electrodes	ES06	Extended 3mm Ball, 5/box	B0312	9006225	909038	60-5181-102	E1563	6156901	DBL-311	31-5000
	ES07	Extended 5mm Ball, 5/box	B0512	9006230	909003	60-5181-103	E1564	6157901	DBL-511	31-5002
	ES20	3/16" Ball, 25/box				138103	E1550			
	ES21	5/32" Ball, 25/box			909136					
	ES50	1/8" Ball, 25/box								
LLETZ Loop Electrodes	ES08	5mm x 5mm Loop, 5/box						6133901	DLP-T11	31-3056
	ES09	10mm x 10mm Loop, 5/box	R1010	9006200	909007		E1559	6134901	DLP-S11	31-5072
	ES10	15mm x 8mm Loop, 5/box	R1507 *	9006205	909011			6135901		
	ES11	15mm x 10mm Loop, 5/box				60-5180-102		6136901	DLP-M11	31-5074
	ES12	20mm x 8mm Loop, 5/box	R2008	9006215	909013			6137901		31-5090
	ES13	20mm x 15mm Loop, 5/box	R2008	9006215	909009			6137901		31-5096
	ES31	20mm x 20mm Loop, 5/box	R2020							31-5098
	ES42	20mm x 12mm Loop, 5/box		9006220 *				E1561	DLP-W11	31-5094
	ES43	15mm x 15mm Loop, 5/box	R1515	9006210						
	ES44	15mm x 5mm Loop, 5/box	R1505		909134					
	ES45	13mm x 13mm Loop, 5/box								
	ES46	10mm x 5mm Loop, 5/box								
	ES47	25mm x 10mm Loop, 5/box	R2510							
	ES48	10mm x 3mm Loop, 5/box								
	ES49	20mm x 10mm Loop, 5/box	R2010		909132	60-5180-104				

			COOPER	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CIRCON	UTAH	MILTEX
LLETZ Loop Electrodes	ES51	13mm x 8mm Loop, 5/box				60-5180-101 *				
	ES52	10mm x 8mm Loop, 5/box			909133*					31-5054
	ES53	15mm x 12mm Loop, 5/box					E1560			
LLETZ Square Electrodes	ES14	5mm x 5 mm Square, 5/box					E1565	6153901		
	ES15	10mm x 4mm Square, 5/box	S1004	9006235			E1566	6154901	DLP-SQ2	
	ES16	10mm x 8mm Square, 5/box			909017		E1567	6155901	DLP-SQ1	
	ES41	5mm x 10mm Square, 5/box								
Arthroscopic Electrodes	AR00	Electrode Hook 90, 5/box					88-1510		E1510	
	AR01	Electrode Hook 45, 5/box					88-1511		E1511	
	AR02	Electrode Acromioplasty, 5/box					88-1512	CS012A	E1512	
	AR03	Electrode Meniscectomy, 5/box					88-15103	CS010A	E1513	

* shows 1mm difference in size

ARTHROSCOPIC ELECTRODES

The Aaron Arthroscopic electrodes are available in four tip styles and designed to fit all standard electro-surgical pencils.



SUPERCUT™ NEEDLES

SuperCut™ Needles (ES 60-63) are made of Tungsten and keep their shape throughout the procedure.

*STERILE: Individually packaged, sterile