

Brewer

135 Series Saddle Seats


Exclusive 360° Dynamic Motion technology optimizes comfort

Designed to provide unmatched comfort

- Innovative 360° Dynamic Motion capability engages core muscles, helping strengthen back muscles over time – reducing back pain and sitting fatigue
- Saddle design eases load on back muscles by tilting pelvis forward to enhance posture, naturally align vertebrae and disks, and relieve shoulder and neck tension
- Enables dentists and hygienists to be more comfortable and productive throughout their workday

Personalized comfort

- Choice of saddle seat sizes and cylinder heights enables you to customize the fit for all users under 5' 3", those that are 5' 3"-5' 11", and those over 5' 11"
- Conveniently adjust cylinder height to exact user preference



360° Dynamic-Motion technology allows the freedom of movement that helps engage core muscles. Adjustable design enables you to limit the range of motion initially. Then, gradually increase the range of motion over time to continue strengthening core muscles, which supports enhanced comfort, daily productivity and career longevity.


Brewer

135 Series Saddle Seat

Brewer exclusives help optimize a practice's bottom-line

Brewer Differentiators

- 360° Dynamic Motion technology promotes back health
- Conveniently adjust the level of dynamic motion desired to suit personal preference
- Standard 5-year warranty



Ergonomically shaped saddle tilts the pelvis naturally to improve spinal posture and enhance long-term comfort. While Brewer's 360° Dynamic Motion technology contributes to personal comfort, daily productivity and career longevity.

135 Amazone Series

All 135 Amazone Seats feature:

- Narrower 13" upholstered seat design
- Press on ring for easy height adjustment
- Turn ring to conveniently adjust level of dynamic motion
- UltraLeather upholstery available in Dove Gray, Raven Wing and Papyrus
- 23" diameter cast aluminum base w/large castors
- 250 lb. weight capacity



135 Jumper Series

All 135 Jumper Seats feature:

- Standard 17" upholstered seat design
- Press on ring for easy height adjustment
- Turn ring to conveniently adjust level of dynamic motion
- UltraLeather upholstery available in Dove Gray, Raven Wing and Papyrus
- 23" diameter cast aluminum base w/large castors
- 250 lb. weight capacity



135 Amazone Models*	Adjustable Height Range	Recommended User Height
■ 135ASG	24"-32"	5' 3"-5' 11"
■ 135ASB	24"-32"	5' 3"-5' 11"
■ 135ASP	24"-32"	5' 3"-5' 11"
■ 135ASG-L	21"-27"	Under 5' 3"
■ 135ASB-L	21"-27"	Under 5' 3"
■ 135ASP-L	21"-27"	Under 5' 3"

135 Jumper Models*	Adjustable Height Range	Recommended User Height
■ 135JSG	24"-32"	5' 3"-5' 11"
■ 135JSB	24"-32"	5' 3"-5' 11"
■ 135JSP	24"-32"	5' 3"-5' 11"
■ 135JSG-H	27"-37"	Over 5' 11"
■ 135JSB-H	27"-37"	Over 5' 11"
■ 135JSP-H	27"-37"	Over 5' 11"

* G - Dove Gray B - Raven Wing P - Papyrus

The Power to Advance. Safely. Financially. Competitively.

www.brewercompany.com | P: 1.888.273.9371

The Brewer Company is an FDA registered & ISO 13485:2003 certified company.

LITD-025-Rev B