



Surgery/ED/ASC/ Trauma Lighting

During the most critical moments, you need a light you can trust.



burton®

Surgery/ED/ASC/ Trauma Lighting

A variety of configurations with easy-to-use features for reliable lighting in critical moments.

From no-drift arms to pure white light with minimal heat output, Burton lights have the features you need for ease of use, versatility and patient comfort. Sleek designs, multiple configuration options and high value on all Burton lights mean you can select the right lights for your facility. LED options provide additional cost and energy savings.



Burton proudly stands behind the reliability of its lighting products with an industry-leading 5-year limited warranty.*

Dept./Specialty	Good	Better	Best
Hospital			
G OR	AIM-200 OR (D) / ALED100 (D)	ALED100 (S)	ALED100 (D)
Trauma	AIM LED (D) / ALED070 (D)	AIM-100 (D) / ALED070 (D)	ALED100 (D)
ED	Outpatient II / CoolSpot II	AIM LED (D) / ALED070 (D)	AIM-100 (S) / ALED100 (S)
ASC			
Endo	AIM LED / AIM-50 / ALED070	AIM-100 / ALED070	ALED100
Plastic	AIM LED / AIM-50 / ALED070	AIM-100 / ALED070	ALED100
Optha	Outpatient LED / Outpatient II	AIM LED / ALED070	AIM-100 / ALED070
Ortho	AIM LED / AIM-50 / ALED070	AIM-100 (D) / ALED070 (D)	AIM-200 OR (D) / ALED100 (D)
G OR	AIM-100 / ALED070	AIM-200 OR / ALED100	ALED100

Examination Lights

AIM® HI LED

- AIM HI 100 - 100,000 lux (9290 fc) at 39.37" (1 m)
- AIM HI 70 - 70,000 lux (6503 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of >95
- 50,000 hours average rated bulb life



AIM-50®

- 45,000 lux (4181 fc) at 39.37" (1 m)
- 3100 K color temperature
- CRI of 96
- Three 50-watt quartz halogen bulbs
- 2,000-hour average rated bulb life



AIM® LED

- 45,000 lux (4181 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of 92
- 50,000-hour average rated bulb life



NOVA EXAM® LED

- 25,000 lux (low level) (2323 fc) at 18" (46 cm)
- 50,000 lux (high level) (4645 fc) at 18" (46 cm)
- 4600 K color temperature
- CRI of 93
- 40,000-hour average rated bulb life



SUPER BRIGHT SPOT®

- 56,000 lux (5185 fc) at 18" (46 cm)
- 3400 K color temperature
- CRI of 98
- 20-watt quartz halogen bulb
- 2,000-hour average rated bulb life



SUPER EXAM® 50 LED

- 39,000 lux (3623 fc) at 24" (61 cm)
- 4000 K color temperature
- CRI of 92
- 40,000-hour average rated bulb life



OUTPATIENT® II

- 86,000 lux (7990 fc) at 24" (61 cm)
- 3300 K color temperature
- CRI of 94
- Three 50-watt halogen bulbs
- 2,000-hour average rated bulb life



OUTPATIENT® II/COOLSPOT® II COMBO

- Outpatient II:
 - 86,000 lux (7990 fc) at 24" (61 cm)
 - 3300 K color temperature
- CoolSpot II:
 - Up to 60,000 lux (5574 fc) at 39.37" (1 m)
 - 4300 K color temperature



Surgical Lights

AIM-100®

- 63,000 lux (5852 fc) at 39.37" (1 m)
- 3500 K color temperature
- CRI of 97
- Three 35-watt quartz halogen bulbs
- 2,000-hour average rated bulb life



AIM-200® OR

- 90,000 lux (8161 fc) at 39.37" (1 m)
- 3700 K color temperature
- CRI of 96
- Larger and adjustable light pattern



APEX

- 120,000 lux (11,148 fc) at 39.37" (1 m)
- 4200 K color temperature
- CRI of 95
- Unlimited arm rotation



Find additional light products at burtonmedical.com

Ordering Information

AIM HI LED	
ALED070FL	ALED070 Floor Stand
ALED070SC	ALED070 Single Ceiling
ALED070DC	ALED070 Dual Ceiling
ALED070W	ALED070 Wall
ALED100FL	ALED100 Floor Stand
ALED100SC	ALED100 Single Ceiling
ALED100DC	ALED100 Dual Ceiling
ALED100W	ALED100 Wall

AIM-50	
A50FL	AIM-50 Floor Stand
A50SC	AIM-50 Single Ceiling Mount
A50DC	AIM-50 Double Ceiling Mount
A50W	AIM-50 Wall Mount

AIM LED	
ALEDFL	AIM LED Floor Stand
ALEDSC	AIM LED Single Ceiling Mount
ALEDDC	AIM LED Double Ceiling Mount
ALEDW	AIM LED Wall Mount

NOVA EXAM LED	
NXFL	Nova Exam LED Floor Stand
NXT	Nova Exam LED Table mount
NXW	Nova Exam LED Wall Mount
NXC	Nova Exam LED Ceiling Mount

OUTPATIENT LED	
OPLEDFL	Outpatient LED Floor Stand
OPLEDSC	Outpatient LED Single Ceiling Mount
OPLEDDC	Outpatient LED Double Ceiling Mount
OPLEDW	Outpatient LED Wall Mount

OUTPATIENT II (COMPLETE MODELS ONLY)	
OP216FL	Outpatient II Floor Stand
OP216SC	Outpatient II Single Ceiling Mount
OP216DC	Outpatient II Double Ceiling Mount
OP216W	Outpatient II Wall Mount
OP216ST	Outpatient II Single Fastrac Mount
OP216DT	Outpatient II Double Fastrac Mount

OUTPATIENT II / COOLSPOT II COMBO	
OC516CC	Combo, Ceiling Mount
OC516CT	Combo, Fastrac Mount

AIM-100	
A100FL	AIM-100 Floor Stand
A100SC	AIM-100 Single Ceiling Mount
A100DC	AIM-100 Double Ceiling Mount
A100W	AIM-100 Wall Mount

AIM-200 OR	
A200FL	AIM-200 Floor Stand
A200SC	AIM-200 Single Ceiling Mount
A200DC	AIM-200 Double Ceiling Mount
A200W	AIM-200 Wall Mount

APEX	
APXSC8	APEX 8' Single Ceiling Mount
APXSC10	APEX 10' Single Ceiling Mount
APXDC8	APEX 8' Double Ceiling Mount
APXDC10	APEX 10' Double Ceiling Mount

NOTE: 120V standard; Additional voltage models also available. Contact Burton for additional information.

* To register your Burton light(s), please visit:
<http://www.burtonmedical.com/warranty-registration/>

Find additional light products and accessories at burtonmedical.com



© 2018 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton[®]

Burton Medical, LLC
 2300 W Windsor Ct., Suite C
 Addison, IL 60101 USA
 Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com