

chattanooga®

THERATHERM® DIGITAL MOIST HEAT PACK

Proven in the Clinic, Portable for Home

The Theratherm® Digital Moist Heat Pack provides moist heat therapy for pain management.¹ Its flannel cover draws moisture from the air and as the packs temperature rises, the moisture is released onto the skin, speeding relief to your patients. Its unique, easy-to-use digital control is convenient and safe.

KEY FEATURES

- Convenient preset time and temperature with on board auto shut-off safety feature and audible buzzer
- New larger easy-read digital hand control for operation, programming, temperature and/or treatment time with new back lit display with one button activation
- Unique and responsibly sourced porous sand retains moist heat while conforming to the body with no chemicals
- Lockout mode protects against inconsistent setting changes throughout treatment sessions
- Moist heat therapy temperatures range from 113° to 167°F (45° to 75°C)
- Treatment sessions up to 60 minutes
- Specially designed flannel cover for efficient moisture retention
- One-year warranty

INDICATIONS

The Theratherm® Digital Moist Heating Pad is intended to provide moist heat therapy for body surfaces.¹ It is capable of maintaining an elevated temperature during use, and temporarily relieving minor muscular or joint pain.¹

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
1031-B	Medium - 13" x 16" (330 mm x 406 mm)
1032-B	Standard - 14" x 27" (356 mm x 686 mm)
1033-B	Shoulder/Neck - 22" x 21" (559 mm x 533 mm)



T 800.336.6569 | F 800.457.4221
1595 Sea Otter Place | Suite 200
Carlsbad, CA 92010 | U.S.A.
djoglobal.com/chattanooga

¹ Theratherm® Digital Moist Heating Pad Instruction Manual (Chattanooga®)

Copyright © 2021 by DJO, LLC - MKT00-11797 Rev A
Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

