


# Clarity Blood Glucose Monitoring System




- ✓ Meter for Single Patient and Multi-Patient Use
- ✓ OTC Approved
- ✓ Strip Ejector
- ✓ Use with Capillary Whole Blood Samples
- ✓ Premeal/Postmeal Indicator



#### Meter Specifications:

- Memory Capacity: 300 Results (Glucose Level, Date & Time)
- Test Specimen: Capillary Whole Blood
- Sample Volume: 0.5 $\mu$ l
- Test Time: 5 Seconds

\* See reverse side for more information.



For more information please email  
[sales@claritydiagnostics.com](mailto:sales@claritydiagnostics.com) or call 877-722-6339  
[www.claritydiagnostics.com](http://www.claritydiagnostics.com)

We Support  
Better Outcomes<sup>™</sup>

# Clarity Blood Glucose Monitoring System



The Clarity Blood Glucose testing system utilize Glucose Oxidase technology that prevents interference with substances such as Xylose, Galactose, and Maltose. Clarity Blood Glucose tests are 15015197 compliant, CLIA Waived, and OTC Approved.

These tests are intended for use in physician's offices, hospitals, at health fairs, and even in the home as an aid in checking blood glucose levels in whole blood to properly manage certain diseases, such as diabetes mellitus.

<p><b>Product Code:</b> <b>CD-BG1</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• BG1000 Glucose Meter ONLY</li> </ul>	<p><b>Product Code:</b> <b>CD-BG4</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• 1 Each of Normal and High Control</li> <li>• Package Insert</li> </ul>	<p><b>Product Code:</b> <b>CD-BG5</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• 50 Test Strips</li> <li>• Package Insert</li> </ul>	<p><b>Product Code:</b> <b>DTG-GL6</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• Lancing Device</li> </ul>
<p><b>Product Code:</b> <b>DTG-GL7</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• Single Use Twist Top Sterile 30G Lancets, 100/BX</li> </ul>	<p><b>Product Code:</b> <b>DTG-SL23100</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• Sterile Single Use 26G Safety Lancets, 100/BX</li> </ul>	<p><b>Product Code:</b> <b>DTG-SL23200</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• Sterile Single Use 26G Safety Lancets, 200/BX</li> </ul>	<p><b>Product Code:</b> <b>CLA-ALSWB100</b></p> <p><b>Contents:</b></p> <ul style="list-style-type: none"> <li>• Sterile Single Use, Medium Size Alcohol Swabs, 100/BX</li> </ul>



For more information please email [sales@claritydiagnostics.com](mailto:sales@claritydiagnostics.com) or call 877-722-6339

[www.claritydiagnostics.com](http://www.claritydiagnostics.com)

**We Support  
Better Outcomes™**