

NE

Clarity Platinum Urine Analyzer

Placement Options and Promotions Available





- Made in the USA Inventory Readily Available
- Correlates 99% to the Siemens CLINITEK Status®+
- Only 510k Approved
 Urine Analyzer Under New
 2018 FDA Guidelines
 - Technology to Prevent Use of Expired Strips
- Significant Cost Savings

- Only Analyzer that Detects Defective Individual Test Pads
- Barcode Scanner Included
- Full Color Touch Screen
- HL7 Compliant

Urine Analyzer with other Urine Analyzers Acon Mission[®] Clarity Platinum

Clinical Comparison of the Clarity Platinum

Clarity Platinum

Parameter	Acon Mission® Analyzer % Agreement with Siemens CLINITEK 500® -Exact Match	Clarity Platinum Urine Analyzer % Agreement with Siemens CLINITEK Status®+ -Exact Match
Glucose	86.3%	98.43%
Bilirubin	92.9%	98.82%
Ketone	88.9%	98.43%
Specific Gravity	66.8%	82.71%
Blood	86.8%	98.04%
Protein	86.3%	97.25%
Urobilinogen	89.5%	99.61%
PH	62.5%	90.57%
Nitrite	98.0%	99.41%
Leukocyte	79.6%	97.84%

For more information call 1-877-722-6339 or email Sales@ClarityDiagnostics.com

- 1. The Acon Mission[®] Analyzer is a registered trademark of Acon Laboratories, INC.
- 2. The Siemens CLINITEK Status®+ as well as the Siemens CLINITEK 500® are registered trademarks of Siemens®.
- Whether the analyzer used for comparison is Siemens CLINITEK Status*+ or Siemens CLINITEK 500*, the comparison was performed using the same Siemens Multistix* 10SG Urine Strips.

We Support Better Outcomes[™]