

Lab X Series, Bariatric, Blood Drawing Chair with Padded Flip Arm and Drawer

MODEL 66002B

- Padded, true, bariatric seat features extra width and depth for added comfort
- Comfortable padded arched backrest
- Height and depth adjustable, padded flip arm
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- A full 35" wide seat (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Extra, steel, center support
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 700 lbs. combined load capacity under normal use (317.8 kg)



MODEL 66002B SPECIFICATIONS

Seat	35" W x 21" D x 20" H
Overall	53" W x 29" D x 36" H
Arm Height Range	27" — 33"
Arm Adjustment	6"

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Standard Upholstery Colors

