





BUDS



CONES



ADVANCED PORTABLE CRYOSURGERY

HISTOFREEZER® FLEX IS THE NEW COLD STANDARD!

Our NEW design works with BUDS and CONES, our colder cryogen will deliver better first-time results, and our returnable canister is the most environmentally friendly canister ever produced for cryosurgery.



ADVANCED PORTABLE CRYOSURGERY

Breakthrough Features Include:



Uses CONES **OR** BUDS



High Efficacy



Single-Handed Design



Effective Freezing



Immediately Cold



Higher Safety




Environmentally Friendly







CLOSED VS. OPEN CONE DESIGN FREEZING TEMPERATURE OVER TIME



Closed cones deliver and maintain a critical freezing temperature (-50C) versus open cones.



Item	Kit Product	Histofreezer® <i>FLEX</i>
No.	Code	Product Description
		ALL KITS CONTAIN:

One Release Clip, One Cone Selection Guide, One Practice Pad, One Instructions for Use and One Pre-paid Return Postage Label

200- 5052	50SM Mixed Kit	Kit Contains: 50 Treatments 1 Canister of cryogen 30 Buds 40 Cones (3, 5, 7, 9, and 12mm) 4 Cone Adapters
200- 6080	80SM Mixed Kit	Kit Contains: 80 Treatments 1 Canister of cryogen 50 Buds 40 Cones (3, 5, 7, 9, and 12mm) 4 Cone Adapters
200- 4051	50SC Cones Only	Kit Contains: 50 Treatments 1 Canister of cryogen 80 Cones (3, 5, 7, 9, and 12mm) 4 Cone Adapters
200- 3050	50SB Buds Only	Kit Contains: 50 Treatments 1 Canister of cryogen 50 Buds
200- 2000	240LR Replacement Kit	Kit Contains: 1 222mL/244g canister of cryogen 1 Sample Pack** Including: 2 Buds/4 Cones/1 Cone Adapter







50 Count Mixed Kit (#200-5052)

Contact your MedPro Account Manager TODAY!