

# Flat Digital Scale

## KEY FEATURES

The *Auto-On* feature of our DS500 switches the scale "on" when the patient taps on the heavy duty glass platform. A patient's accurate weight measurement is easy-to-read because of the unit's 2" high back-lit LCD display. The durable tempered glass weighing surface is 12" x 12" and the scale's sleek contemporary design is just 3/4" high. The DS500 can display weight in lbs or kg.

- Auto-On function for easy start up
- Spacious weighing platform
- Convenient 0.75" scale height
- Durable tempered glass construction
- Easy-to-read LCD display with 2" numerals
- Quickly switch between lb or kg

## SPECIFICATIONS

<b>Warranty:</b>	Two Years
<b>Power Input:</b>	3V Lithium Battery (Included)
<b>Battery Life:</b>	3,000 Weighments
<b>Maximum Capacity:</b>	400 x 0.2 lb (180kg x 0.1kg)
<b>Accuracy:</b>	+/- 0.2 lb ( 0.1 kg)
<b>Construction:</b>	Heavy Duty Tempered Glass
<b>Dimensions:</b>	11.75" L x 12.25" W x 0.75" H
<b>Displayed Units:</b>	lb, kg

