



WE CARE LIKE FAMILY™

Black Arrow™ Examination Gloves

Make Your Mark

Dynarex Black Arrow™ Latex Examination Gloves are ideal for law enforcement, correctional facilities, tattoo shops, and high-end beauty salons. They are also a good choice for light-duty maintenance when protection is required.

- *Maximum protection, puncture resistant, durable, tactile sensitivity, and comfortable fit*
- *Powder free to minimize skin irritation (≤ 2 mg residual powder per glove)*
- *Textured for secure grip in wet or dry use*
- *Polymer coated and ambidextrous*
- *3-5 mil. range min. thickness at finger & palm*

Technical Specifications: Dynarex Black Arrow™ Latex Examination Gloves are medical grade gloves that meet or exceed ASTM D3578 specifications. Protein content is 50µg/dm² or less.

CAUTION: This product contains natural rubber latex which may cause allergic reactions. Do not expose this product to any person known or suspected to be latex sensitive. Do not use in latex-free areas or on patients with spina bifida.



Reorder No.	Black Arrow™ – Powder Free	Quantity
2321	Black Latex Exam Gloves - S	10/100/Cs
2322	Black Latex Exam Gloves - M	10/100/Cs
2323	Black Latex Exam Gloves - L	10/100/Cs
2324	Black Latex Exam Gloves - XL	10/100/Cs

See our website for our Other Types of Gloves



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com