

LONGTERM CARE BEDS

D300 Series Low Beds



WE CARE LIKE FAMILY™

FEATURES AND BENEFITS

- Constructed from welded steel
- Features true vertical lift
- Deck to floor range 7.8"—30"
- Ultra-quiet DC motor has 5 functions: head, foot, bed height, Trendelenberg, and Reverse Trendelenberg adjustment
- Integrated, tool-free, length extension extends: 80"—84"—88"
- Equipped with eight 3" casters
- Wall bumper guard protection
- Weight Capacity 500 lb.



Warranty: 15 years on structural steel, 3 years on electrical & mechanical, 1 year on all other parts and components from original date of purchase.

COMPARE TO:

Prime Care P703 by Drive DeVilBiss Healthcare

ALSO AVAILABLE:

- Tool-free, Removable Headboards/Footboards:
Available in Mahogany, Amber Cherry, and Light Oak or Composite Plastic
- Metal Swing Rail, Countoured Swing Rail, or Assist Bar
- Footboard with Embedded Nurses Control
- Battery Backup
- Extensions to Expand D300 Bed and Headboard to 42"

Nurses Control Panel



Series	Height		Expansion		Extended Heel Lift	Central Locking System	Weight Capacity
	Min	Max	Length	Width			
D300	7.8"	30"	80"—84"—88"	36"—42" *	✓	✓	500 lb.

* with optional headboard



Find Us On:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR
Call us at: (888) 396-2739 • www.dynarex.com