




Elastic Bandage with Self Closure



Self Closure Feature



Ideal for Low Compression Situations

These premium elastic bandages with self closure have a 1" self closure start strip and a 2" self closure strip for easy and secure fastening. Premium Elastic Bandages are a non-slip compression and support wrap that uses no pins, clips, or tape. For best results, stretch and apply the elastic bandage before securing the self closures.

- Latex Free
- These bandages are tannish in color (white with tan)
- Easy-to-use, fast and apply maximum compression
- Comfortable and strong enough for laundering and repeated usage
- Premium elastic bandages are convenient and efficient with the use of self closure feature

Reorder No.	Description	Quantity
3658	2" Elastic Bandage	5/10/Cs
3659	3" Elastic Bandage	5/10/Cs
3660	4" Elastic Bandage	5/10/Cs
3661	6" Elastic Bandage	5/10/Cs



Visit our website to see our ActiSplints (3353, 3528-3530).



Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.

For more information contact your local Dynarex® distributor sales representative. Or call 888-396-2739 to find your nearest distributor.

Already mobile? Scan your smart phone here for more information.

