



WE CARE LIKE FAMILY



## Laryngoscope Blades /Handles

**Used as a primary tool for examining the interior of the larynx or for placement of an endotracheal tube**

- Variety of blade attachments (Miller and Macintosh) with fiber optic light source
- Eliminates the risk of cross-contamination from both blades and handles
- Stainless steel construction
- Laryngoscope blades are rigid
- Precise control to elevate the epiglottis
- Tip has non-pinching design
- Blades are sterile; handles are non-sterile
- Blades fit on both Penlight and Standard handles
- Laryngoscope blades and handles are interchangeable with Welch Allyn



4024-4028  
Macintosh Blade



4029-4034  
Miller Blade




4035  
Disposable Handle  
(4036 not pictured)

Reorder No.	Description	Quantity
4024	Macintosh Blade #0 – Disposable, Fiber Optic	10/Pkg
4025	Macintosh Blade #1 – Disposable, Fiber Optic	10/Pkg
4026	Macintosh Blade #2 – Disposable, Fiber Optic	10/Pkg
4027	Macintosh Blade #3 – Disposable, Fiber Optic	10/Pkg
4028	Macintosh Blade #4 – Disposable, Fiber Optic	10/Pkg
4029	Miller Blade #00 – Disposable, Fiber Optic	10/Pkg
4030	Miller Blade #0 – Disposable, Fiber Optic	10/Pkg
4031	Miller Blade #1 – Disposable, Fiber Optic	10/Pkg
4032	Miller Blade #2 – Disposable, Fiber Optic	10/Pkg
4033	Miller Blade #3 – Disposable, Fiber Optic	10/Pkg
4034	Miller Blade #4 – Disposable, Fiber Optic	10/Pkg
4035	Disposable Handle – Standard*	10/Pkg
4036	Disposable Handle – Penlight**	10/Pkg

\*Standard handle uses (2) "C" batteries.  
\*\*Penlight handle uses (2) "AA" batteries.  
Batteries not included.

Welch Allyn is a registered trademark of Welch Allyn®

See our website for our Endotracheal Tube Introducers (Bougie) (#4582-4583)



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)