



**WE CARE LIKE FAMILY**

## Medicine Cups

**Manufactured of durable material with easy-to-read gradations on all sides**

### Item 4252

- Perfect for dispensing both liquids and medications
- Constructed of virgin polypolyplene resin - with no lip
- Calibrated for measurements: cc or ML: (2.5, 5, 7.5, 10, 15, 20, 25, 30)
- Can hold up to 1 fluid ounce

### Item 4258 (Economy)

- Calibrated for measurements: cc or ML: (2.5, 5, 7.5, 10, 15, 20, 25, 30)
- Constructed of virgin polypropylene - with a lip
- Can hold up to 1 fluid ounce



4252



4258

Reorder No.	Description	Quantity
4252	Medicine Cup-1 oz.	50/100/Cs
4258	Medicine Cup-1 oz. (Economy)	50/100/Cs

See our website for our Disposable 3 oz. Cups (#4235)



Find us on:



**dynarex**® is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)