

COMPRESSOR NEBULIZERS

Dynarex Resp-O₂™ Compressor Nebulizers feature compact, ergonomic designs that allow for the easy inhalation of prescribed medication. These nebulizers help users to efficiently manage respiratory ailments.

- Small Volume Nebulizer (SVN)
- Converts prescribed medication into an aerosol mist for easy inhalation
- Quiet operation levels
- Includes reusable nebulizer cup
- Not made with natural rubber latex or DEHP

INPUT	AC 120 V, 60 HZ
EXTREME PRESSURE	>205 kpa
POWER CONSUMPTION	≤200 VA
MEDICATION CAPACITY	10 mL
PARTICLE SIZE	≤5 μm & 60%
MMAD	2.0 μm
AVERAGE NEBULIZATION RATE	≥0.2 mL/min
NOISE LEVEL	≤ 60 dB
COMPRESSOR PRESSURE RANGE	30-58 psi
OPERATION PRESSURE RANGE	14.5-27 psi (1 bar-1.9 bar)
LITER FLOW RANGE	16-23.5 lpm
ACCESSORIES	Adult Mask, Pediatric Mask, Reusable Nebulizer Cup, Mouthpiece, T-Piece, 7' (2.1 m) Tube, 5 Filters



Elite

REORDER NO. 34400

DIMENSIONS
 29 × 18 × 11 cm

PACKAGE QUANTITY
 4/Case



Compact

REORDER NO. 34401

DIMENSIONS
 15.5 × 13 × 8.3 cm

PACKAGE QUANTITY
 6/Case



Compact (International)

REORDER NO. 34408

DIMENSIONS
 15.5 × 13 × 8.3 cm

PACKAGE QUANTITY
 6/Case

